




1 x 3 30ml syringe, 20ml samples, 100g batches




1 x 4 30ml syringe, 24ml samples, 100g batches




1 x 5 30ml syringe, 20ml samples, 100g batches




1 x 6 30ml syringe, 20ml samples, 100g batches




1 x 7 30ml syringe, 24ml samples, 100g batches




1 x 8 30ml syringe, 21ml samples, 100g batches




1 x 9 30ml syringe, 24ml samples, 100g batches




1 x 10 30ml syringe, 27ml samples, 100g batches




1 x 11 30ml syringe, 20ml samples, 100g batches




3 x 3 30ml syringe, 20ml samples, 100g batches




3 x 5 30ml syringe, 24ml samples, 200g batches




3 x 6 30ml syringe, 24ml samples, 200g batches




4 x 3 30ml syringe, 24ml samples, 100g batches




4 x 4 30ml syringe, 27ml samples, 200g batches




5 x 4 60ml syringe, 48ml samples, 300g batches




5 x 5 60ml syringe, 48ml samples, 300g batches



6 x 4 60ml syringe, 45ml samples, 300g batches



6 x 5 60ml syringe, 40ml samples, 300g batches



7 x 5 60ml syringe, 48ml samples, 300g batches

